

# Technical documentation

---

Product name:	<b>qRE <i>Rhodiola rosea</i> L., roots</b>
Substance:	<i>Rhodiola rosea</i> L., roots dry extract
Plant source common names:	en: Roseroot; fr: Rhodiole
Reference:	E0059
Packaging:	100 mg in a 1.5 ml borosilicate amber vial
Storage conditions:	Keep container closed. Protect from light and moisture. Keep inferior to -15 °C.
Retest:	12 months

## Botanical identification of plant source

---

Plants in our botanical garden are identified and a herbal voucher is prepared by an expert botanist.  
Each batch collected for extraction is verified and identified.

**Reference:** Flora Europaea, Cambridge, University Press, 1993, Vol 1, p 363

## Method of production of dry extract

---

Whole plant or plant parts are collected, freeze-dried and coarsely ground. Extraction is performed by maceration in 50 % (v/v) aqueous ethanol for 48 hours at room temperature. Ethanol is then evaporated under reduced pressure at less than 40 °C and the aqueous residue is freeze-dried.

## Organoleptic characteristics of dry extract

---

Colour: Light brick red

Odour: Characteristic

Form: Fine powder

## Recommended methods for use

---

Weight a precise weight of qRE and solubilise in the recommended solvent at the concentration indicated in the HPLC or HPTLC method described in this document.

Sonicate for 90 seconds (70 W). Filter on a 0.45 µm PVDF membrane and put the resulting solution into HPLC dispenser or apply on the HPTLC plate.

Dose and analyse your extract with qRExtract using the HPLC / HPTLC methods described in this document or using your own methods.

## HPTLC

### Detection of rosavin

**Layer:** 10 × 10 cm HPTLC Nano-Sil-20 UV 254 (Carl Roth ref. N084.1)

**Thin layer conditionnement:** 1 h at room temperature and 33 % relative humidity

<b>Elution solvent:</b>	<u>Elution solvent compound</u>	<u>Volume (ml)</u>
	ethyl acetate	77
	methanol	15
	H <sub>2</sub> O	8

**Developing distance:** 70 mm from the lower edge


**Initial spot volume and concentration:**

qRE:	5 µl of a 2.6% (w/v) solution in 50 % aqueous ethanol
rosavin:	2 µl of a 0.1% (w/v) solution in 50 % aqueous methanol

**Reagent mixture:**

Aniline - diphenylamine-phosphoric acid reagent

Preparation: Dissolve 1 g of diphenylamine in 40 mL of acetone, add 1 mL of aniline, and mix. Carefully add 7.5 mL of phosphoric acid, and mix. Dip in the reagent mixture and dry for 5 minutes at 120 °C. Expose to visible light.



## HPTLC

### Detection of salidroside

**Layer:** 10 × 10 cm HPTLC Nano-Sil-20 UV 254 (Carl Roth ref. N084.1)

**Thin layer conditionnement:** 1 h at room temperature and 33 % relative humidity

<b>Elution solvent:</b>	<u>Elution solvent compound</u>	<u>Volume (ml)</u>
	ethyl acetate	77
	methanol	15
	H <sub>2</sub> O	8

**Developing distance:** 70 mm from the lower edge

**Initial spot volume and concentration:**

qRE: 1,5 µl of a 2.6 % (w/v) solution in 50 % (v/v) aqueous ethanol

salidroside: 1 µl of a 0.1 % (w/v) solution in 50 % (v/v) aqueous methanol

**Reagent mixture:** Iron(III) chloride - Potassium ferricyanide


Preparation:

1) dissolve 4.5 g of ferric chloride in 100 mL of 50 % (v/v) aqueous ethanol.

2) dissolve 1g of potassium ferricyanide in 100 mL of 50 % (v/v) aqueous ethanol.

Dip in a combination (1:1) of two reagent mixtures and dry for 10 minutes at room temperature.





Expose to visible light.



## HPLC

<b>Precolumn:</b>	Ascentis® Express C18 0.5 cm × 3.0 mm 2.7 µm		
<b>Column:</b>	Ascentis® Express C18 15 cm × 3.0 mm 2.7 µm		
<b>Sample:</b>	10 µl 1.5 % qRE (w/v) solution in 25 % (v/v) aqueous ethanol		
<b>Flow:</b>	0.45 ml/min		
<b>Temperature:</b>	25 °C		
<b>Mobile phase:</b>	A: 0.1 % formic acid (v/v) in water B: 0.1 % formic acid (v/v) in acetonitrile		
<b>Detection:</b>	Diode Array Detector, 256 nm		
<b>Gradient:</b>	Time (mn)	A %	B %
	0	97	3
	65	79	21
	90	60	40

## Quantified substances

Compound	CAS No	2D Structure	Peak No
Salidoside	10338-51-9		1
Rosavin	84954-92-7		4
Cinnamyl alcohol	104-54-1		5
Rhodosin	86831-54-1		6
Unknown	NA	NA	2, 3, 7, 8